tifix[®] Ilizarov Ringfixateur & Hexapod System

tifix[®] Ilizarov ring fixator & hexapod systems




Prinzip des Segmenttransports beim **Ilizarov**-Ringfixateur
Principles of segment transport using the **Ilizarov** ring fixing device

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Segment transport utilizes the experiences gained in lengthening extremities. Various investigators have proven that continuous or pulsed extension stimulus between two bone ends generates granulation tissues in which fibroblast-like cells form collagen fibrils around central capillaries. Osteoblasts which form osteoid tissues arrange themselves evenly in this soft tissue structure. This results in the formation of a typical long hollow bone. This type of bone formation is termed “intramembranous bone formation” [Pesch and Wagner 1974, Monticelli 1981, Paley 1988, Ilizarov 1989, Aronson et al 1989, Delloye et al 1990, De Pablos and Canadell 1990, Lascombes et al 1991, Shearer et al 1992].

Contrary to extension, segment transport involves cutting a vital hollow bone segment from the vital bone tube and gradually shifting it into the defect that must be filled, using a so-called drawing mechanism. Finally the shifted segment “docks” on the original bone tube and is here immersed under compression or an autogenic spongiosa restoration is carried out in the “docking site” for more rapid intergrowth. At the same time, the transport mechanism is changed into a fixation ring (Fig. 1a – e).

Segment transport can be done in various ways. The easiest method is simple transport (Fig. 2a), which can be placed from the proximal to the distal direction or conversely from the distal to the proximal direction. Since the proximal circulation is more favorable than the distal circulation, this transport direction is favored, wherein the separation of the bone tube which Ilizarov refers to as “corticotomy” should lie in the metaphyseal segment if possible.

In double transport (Fig. 2b), 2 segments are transported against each other, wherein transport fundamentally takes place 1 x from the proximal to the distal direction and 1 x from the distal to the proximal direction (possibly here at half the speed).

The significant advantage of the double transport lies in shortening the required distraction time (= transport time). Tandem transport is another variant; in this, both segments are moved from the proximal to the distal direction or the reverse; however, the more distant segment must be drawn at double the speed (Fig. 2c). Double and tandem transport are extraordinarily susceptible to problems and can only be used in very differentiated patients.

Prinzip des Segmenttransports beim Ilizarov-Ringfixateur
Principles of segment transport using the Ilizarov ring fixing device

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


When using the ring fixing device, this means that every fragment must be stabilized with 2 rings each (Fig. 1 a-d, 3a-d).

Biomechanical aspects are significant for the stiffness of the system:

In the ring pair, it is particularly necessary to note the following:

- a) Ring size (as small as possible)
- b) Number of wires (at least two between 30° and 60° to each other)
- c) Wire tension (500 – 700 N) (the wire begins to slip at 510 N)
- d) Inner distance between the rings (as large as possible, minimum distance 4 cm between the rings)

2. Bone separation which is gentle on tissues and maintains vascularity = corticotomy

[Allan 1948, Anderson 1952, Kawamura et al 1968, Ilizarov 1989, Brutscher et al 1992, Kenwright und White 1993]:

It was of decisive importance for Ilisarow to maintain the medullary vessels. Therefore he and his co-workers developed a special bone separation method wherein the corticalis is ventrally, medially and laterally separated with a special chisel, while the dorsal corticalis was fractured by rotation of the fragments against each other.

In the osteotomy of the bone tube with a chisel or saw, the medullary blood vessels are cut; as Brutscher was able to prove, this does not prevent formation of the intramembranous bone growth, but delays it by approx. 4 weeks.

Prinzip des Segmenttransports beim **Ilizarov**-Ringfixateur
Principles of segment transport using the **Ilizarov** ring fixing device

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After corticotomy and moderate separation of the fragments (max. 2 mm), it is necessary to wait a few days until stretchable granulation tissues have formed between the fragments. Regarding the duration of the distraction phase, recommendations vary from 4 – 12 days. Fundamentally, the following applies:

The more severe the soft tissue and/or bone damage, the longer it is necessary to wait before starting distraction.

4. Extent and rhythm of the distraction

[Ilizarov 1989, White and Kenwright 1990]:

Ilizarow and his co-workers in particular investigated the distraction rhythm and the distraction distance in extensive animal tests. Herein the particularly successful circumstances were shown to be those which had not been used in the decades-old extension techniques of the Western world, due to which these techniques had so often shown poor results.

The smaller the distraction step and the more frequently it was done, the better was the resultant bone formation.

In other words, this means that the distraction distance of 1 mm per day which is recognized as particularly effective should be reached in as many individual steps

as possible. For this reason, great significance is given to the development of distraction motors with a very slow drawing speed. The rule of thumb is as follows:

4 x 0,25 mm/day. In children and on the thigh, the distraction distance can be increased to 2 mm/day in individual cases.

5. Maintaining and/or restoring the function of the extremities

[Ilizarov 1992, Schmidt and Wolter 1995]:

Ilisarow recognized and always emphasized the significant importance of functional therapy.

However, many authors neglect this perspective in particular.

It must be emphasized that every success of a segment transport depends on the use therapy which must be applied simultaneously in order to regain a functional extremity.

Durchführung der **Ilizarov**-Ringfixateurmontage Mounting the **Ilizarov** ring fixing device

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Mounting the Ring Fixator on the lower leg

We will explain the mounting form here as we have seen it in its most favorable form after 18 years of constant use.

We now only rarely use assemblies which were produced before surgery because – with some practice – intraoperative assembly as it is shown here can be quickly performed.

We also do not prefer the original Russian mounting form in which all wires are first laid and the ring assembly is then done.

These assemblies are often not very symmetric; this increases the hindrance to the patient.

We are not aware of any serious disadvantages of our assembly.

To mount the Ring Fixator on the lower leg, the patient lies prone; the lower leg is washed up to approximately mid thigh, including the entire foot, and placed in a sterile position.

Herein the thigh and foot are cushioned so that the lower leg itself lies approx. 10 cm over the table.

After disinfection, the so-called landmarks are marked with a water resistant pen in order to make assembly easier and simpler.

One marks the following anatomical structures:

- patella
- knee joint gap, medial and lateral
- fibula protuberance
- tibial tuberosity
- ventral edge of the tibia (rotation control)
- outer ankle bone
- inner ankle bone
- talocalcaneal joint gap
- gap of the large toe base joint
- position of the bone defect

Then the required ring size is selected by holding a ring that has been screwed together over the proximal tibia. The ring is asymmetrically placed.

There must be as little free space as possible medially and ventrally between the inner edge of the ring and the soft tissues.

However this must be at least one finger width - and at least two finger widths laterally and dorsally – in order to avoid cutting off the soft tissues with postoperative swelling.

Caution: Due to the possible swelling, the ring should be chosen so that it is more likely to be too large than too small (biomechanically, a smaller ring is more stable).

Durchführung der **Ilizarov**-Ringfixateurmontage Mounting the **Ilizarov** ring fixing device

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Foot mounting is not required for a simple lower leg mount – with the lower leg defect approximately in the center of the shaft – consequently 4 rings which are approximately equal in size are required for such assemblies (otherwise there are problems with the symmetry of the assembly, or one has to use conical disks) (Fig. 3a-d).

The 2 ring groups are held over the lower leg, the rings are set at right angles to the tibial axis and then the distances between the rings are set (Fig 3a-d) – as already mentioned above – as far as possible apart because this results in the greatest biomechanical rigidity. Then both pairs are connected to the threaded rods which protrude the least (with a ring distance of 7 cm, threaded rods with a length of 10 cm are sufficient). The assembly connections are always made only ventrally and dorsally, wherein the threaded rods are placed directly into the first hole next to the ring lock (Fig. 3c) (this is important because the two threaded rods medially and laterally would otherwise be in the way during assembly when placing the wires, and furthermore the required ring holes might be occupied) (Fig. 3c).

After mounting the two pairs of rings, the length of the distance between the inner rings is measured, and then the connection of both ring pairs is again performed with the shortest possible threaded rods (only select longer rods if an extension is also planned parallel to the segment transport).

Then select the first hole next to the lock on the rings, or opposite the already occupied hole on the opposite side (through ring pair) (Fig. 3c+3d).

This results in a four-ring construction, wherein the two fixing rings are held by two threaded rods and the defect is also bridged by two threaded rods, and all threaded rods are placed directly adjacent to the rods (Fig. 3c+3d).

Now the construction is again held over the lower leg (distance to the joint lines approx.

1.5 – 2 cm in the knee, 1 cm in the upper talocalcaneal joint); the correct position of the defect is taken into account (if a sequestrectomy is also required, the expected defect must be used as the basis) and the distance to the soft tissues is verified. When everything is correct, all nuts are tightened, measuring equal distances ventrally and dorsally, in order to obtain a stable construction.

Durchführung der **Ilizarov**-Ringfixateurmontage
Mounting the **Ilizarov** ring fixing device

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
Then the construction is fastened on the lower leg. For final orientation of the system, taking rotation into account (important!), so-called reference wires without olive are applied cross-wise (which the Americans refer to as a lateral wire) on the respective outer rings – that is, close to the knee and talocalcanean joint – and screwed on to the rings with slotted wire fixation bolts under tension.

Then the position of the ring system on the lower leg is checked (joints unrestricted? rotation correct? movement in the knee and talocalcanean joint unhindered? defect correctly positioned?) (Fig. 3d).

To insert the fixation wires:

Since approx. 1994, we have only used titanium wires with a thickness of 2 mm because their infection-free standing time has been shown to be the longest. The wires have a specifically ground tip (Fig. 4a), so that they can be drilled through the bone with relative ease and little heat generation. When they are inserted through the soft tissues to the bone, the wires are manually pushed – not drilled – and then drilled through the bone (in the case of hard bone, frequent drilling pauses are necessary, otherwise ring sequesters are formed – heat during rapid drilling is up to 600°C -) and after the bone has been traversed, the wires are driven further forward with a hammer, again not drilled (so that no soft tissue structures are wrapped up!).

Fixation takes place with so-called slotted wire fixation bolts. The wire is clamped between the ring and the screw (Fig. 4b, 4d). We pre-tension the wires with tensioning grips which generate between 700 and 1000 N of tension (Fig. 4d). When the wire has been tensioned, the slit screw must be screwed in as far as possible without rotating the head against the wire; if it is too loose, the wire will slip. Unstable assembly is the result. This always causes problems and advances pin infection. If the counter-nut is tightened to 10 Nm, the wire already begins to slip at wire tensions of 510 N; if the counter-nut is tightened to 20 Nm (very firm tightening!), it only begins at 740 N [own biomechanical measurements].

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Mounting the **Ilizarov** ring fixing device

If the ring mount is correctly placed on the lower leg, it is now possible to start introducing the actual fixation wires. We always start proximally.

For extension, correction and segment transport on the lower leg, we only use olive or stopper wires.

The use of these stopper wires is rejected by some authors [Green 1991, Calhoun et al 1992, Catagni (pers. communication)], wherein we cannot share their doubts (more rapid infection, rapid loosening, more complaints). For us, the olive K-wires have proven themselves extremely well in more than 700 mounts.

The olive K-wires guarantee correct position of the tibia in the long term even with short fragments, osteoporotic bone and shaft defects.

The use of bone screws [Green 1991] does not represent a solution for us because the insertion of screws completely alters the biomechanics of the fixing device model.

The success of the ring fixing device is based on the specific rigidity of the system with tensioned wires: the system is relatively soft in the axial direction, however its flexion and torsion rigidity are relatively high.

On the outer rings, 2 olive K-wires per ring are inserted in addition to the reference wire. The one from the lateral direction with an olive lies directly ahead of the fibular protuberance.

The angle to the right lower leg axis is approximately 250°. The degree divisions refer to the right extremity which lies horizontally on the table:

- ventral 0°
- medial 90°
- dorsal 180°
- lateral 270°
- left: angle mirrored)
- the one from the medial direction with olive from the medial direction with approx. 120°.

The tibial wire which comes from the medial direction is significant because this could tangent and therefore block tendons of the pes anserinus.

The wire must be positioned behind the tendons toward the back of the knee. Naturally it must not injure the blood vessels during insertion.

After pushing through the soft tissues, one feels for the medial tibia edge with the tip of the wire, and penetrates the bone here.

Durchführung der **Ilizarov**-Ringfixateurmontage Mounting the **Ilizarov** ring fixing device

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After pulling the wire through to the lateral side, wherein one must not strike the tibial tuberosity, one verifies the unrestricted position of the wire by bending the lower leg in the knee joint. If the wire is correctly positioned, it does not move during the motion. If it tangents a tendon, it jumps during the motion and must be newly placed.

The wire may only be fixed and tensioned after this test (Fig. 5a).

The second ring from the proximal side is occupied with an olive K-wire from the lateral side (with approx. 270°) = lateral wire. (Fig. 5b)

Due to the N. peroneus in the second and third fifths of the lower leg, a fibular wire must not be used here (with approx. 240 – 260°).

From the medial side, an olive K-wire is inserted approximately parallel to the medial facet – therefore the Americans refer to it as the medial facet wire – at approx. 110 - 130° (Fig. 5b).

When pushing through the soft tissues – before drilling through the bone – the crossing muscles must be pre-stretched in the case of all wires which cross muscles; that is, the foot must be kept in maximal dorsal flexion when crossing the M. soleus or M. gastrocnemius medially so that foot function is not later blocked. The assistant must ensure this.

If there is still a wrong position in one or both segments before the ring fixing device is mounted, this malposition is corrected during assembly.

A rotational malposition is corrected when inserting the reference wires.

Varus, valgus, ante- or recurvation malpositions and/or combinations are gradually corrected while mounting the individual rings.

Durchführung der **Ilizarov**-Ringfixateurmontage Mounting the **Ilizarov** ring fixing device

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

Correction is easily done after approximate assembly of the main segments proximally and distally performed by fastening the outer rings proximally and distally as shown above and then performing the correction when mounting the respective inner rings.

For instance, in the case of a varus malposition, a lateral or nearly lateral olive K-wire is brought in on the second ring from the medial side and the wire is gradually tightened with the tensioning grips until the malposition is corrected (Fig. 6a+6b).

Caution: If the wire is to remain in place, it must lastly be tensioned medially in the lateral direction after fixing the slotted fixation bolt. If it disturbs in this position, it can also be removed again after introducing the 2 typical fixation wires.

The distal fourth ring (from the proximal direction) is mounted with 2 olive K-wires next to the reference wire (laterally), wherein the medial wire is placed at approx. 120° (medial facet wire). Herein one seeks the medio-dorsal edge of the pylon with the tip of the wire (caution: avoid slipping into the A. tibialis posterior in the dorsal direction).

One then drills in the direction of 1 to 2 finger widths laterally from the tibialis anterior tendon in order to avoid injuring the A. dorsalis pedis or A. tibialis anterior here.

When drilling these wires, the foot must be held so that the crossing musculature is pre-stretched against the direction of function – important: loss of function! (Fig. 5c+5d)

After this wire has been pre-tensioned and fixed, the third wire in this ring is inserted at approx. 240° through the fibula (fibular wire) (Fig. 5c, 5d).

This prevents the talocalcaneal joint fork from spreading and any gaping of the syndesmosis gap is corrected. The wire aims at the medial surface of the tibia; it must not be inserted so steeply that it laterally passes the tibia or strikes the tendon of the tibialis anterior. When inserting 3 wires on a ring, it is certainly possible to strike other wires with the second or third wire. Further drilling is not possible, or the wire is more easily or strongly steered off its course.

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
If the deviation is greater than 1 cm, the wire should be reinserted. Up to 1 cm, the gap is equalized using long slotted fixation bolts, inserting 1 or 2 thick spacer disks without a slit. Greater distances to the ring must never be ignored and equalized by pushing the wire against the ring. This produces a contour effect.

If it is not possible to insert the wire close to the ring because the other wires are struck, one chooses an entry point which is primarily 1 cm away from the ring, and then fastens the wire with the above mentioned thick spacer disks. Biomechanically, this also results in slightly improved stability!

The third ring is occupied with 2 olive K-wires – the medial wire again at approx. 120° (medial facet wire), and the lateral wire at 240 – 250° (fibular wire), wherein the fibula can - but does not have to - be captured again. In malpositions of the distal fragment/segment, corrections can be made when inserting the wires.

In the end, a stable assembly with a correct axis must result. If this is not the case, further corrections can still be made. For this purpose, the threaded rods between ring 2 (proximal segment) and ring 3 (distal segment) are removed, the correction is made e.g. under image converter control, and the threaded rods are fixed on the rings with so-called conical disks. The conical disks (Fig. 7a+7b) permit fixation in an angled position between the threaded rod and the ring (approx. 15° of play in all directions), while simple fixation between the nut on the threaded rod and the ring always results in a right angled position. This is of significant importance. Even a minor angled position is corrected into a right angled position without conical disks, that is, corrections are not possible without conical disks!

When the ring mounts have been completed, another 2 threaded rods, respectively, are mounted between the rings; they are as evenly and/or harmonically distributed between the two which are already in place as possible.

Due to the wire fastenings, a completely harmonic, even distribution is rarely possible.

Biomechanically, three evenly distributed threaded rods are just as stable as four almost evenly distributed threaded rods.

Durchführung der **Ilizarov**-Ringfixateurmontage Mounting the **Ilizarov** ring fixing device

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Assembly with a short distal segment/fragment or simultaneously existing foot malposition:

If the distal fragment is shorter than 4 cm, a 2 ring mount is not possible because the inner distance of the rings would be noticeably less than 4 cm, which is the minimum distance in biomechanical terms for sufficient stability.

In this case – or if a foot malposition must be corrected at the same time (often contractile pointed supination malposition) the rear foot is included into the lower leg assembly.

The decisive olive K-wires of the rear foot are two crossed, opposing olive K-wires in the heel bone. Herein the medial wire is placed at approx. 110° (fully situated in the body of the heel bone), and the lateral wire is placed at approx. 250° (also fully situated in the body of the heel bone) (Fig. 8a). Seen from the front or back, the wires are shifted by 0.5 - 1 cm so that the ring can be placed; the wires are fixed on this ring from the top or bottom after pre-tensioning (max. 600 N). (Fig. 8a-c)

If – as we prefer – a 5/8 ring is used – which should have approximately the same size as the lower leg rings – this ring is dorsally oriented approximately to the last lower leg ring, so that, insofar as possible, suitable threaded rods can be inserted at right angles between the lower leg assembly and the foot assembly. Herein two threaded rods – dorsomedial and dorsolateral – are used to orient the foot. (Fig. 8b+8c)

If the foot is in the neutral position, two further threaded rods can then be inserted farther in the ventral direction. This results in an extraordinarily stable construction between the distal lower leg and foot which replaces the second ring in the pylon very well with a short fragment and only one ring.

Durchführung der **Ilizarov**-Ringfixateurmontage Mounting the **Ilizarov** ring fixing device

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If the foot is in the neutral position to the lower leg after all four threaded rods are tightened, the upper and lower talocalcaneal joint must be slightly pre-stretched by first detaching the dorsally positioned two threaded rods on one side and then performing a distraction of approx. 0.5 cm (5 turns of the nut).

Then – after fixing the dorsal threaded rods – the ventral rods are correspondingly corrected. This results in a slight dehiscence of the joint surfaces in the upper and lower talocalcaneal joints and prevents ankylosis of the joints in the long term. At the same time, the stability of the joint connections is increased.

The foot is frequently not in a neutral position to the lower leg, but has complex malpositions – usually supination and pointed foot positions – wherein the heel has a varus position.

The heel bone wires are always inserted parallel to the sole of the foot, entirely regardless of how large the angled position of the foot ring to the lower leg is. This is of decisive importance; otherwise correction is not possible. The ring is fastened so that the malposition is imitated. Simple connections between the lower leg assembly and the foot are therefore not possible. (Fig. 8b, 8c)

Slight malpositions are equalized with conical disks; greater malpositions are equalized by including joints in the direction of the correction.

When assembling the threaded end pieces (“female end pieces”), care must be taken to ensure that the threaded rod is screwed to the end of the part receiving the threaded section (Fig. 9a) and not, as shown in Fig. 9b, to the stop, for then it would no longer be possible to mount two end pieces. (Fig. 9c)

Before completing the lower leg/foot assembly with foot malposition, the upper and lower talocalcaneal joints should always be stretched after introducing the four threaded rods or 4 longitudinal bearers with joints, wherein the medial longitudinal bearers must be distracted more strongly than the lateral ones due to varus malposition of the rear foot.

This always takes place step by step; otherwise no correction is generated – only the ring will be angled! With a solid contracture, the correction should take place more slowly. With a solid contracture, stretching is permitted to a maximum of 2 mm – otherwise skin necrosis and/or enormous pain may occur!

With definite malpositions and/or solid contractures, the corrections must absolutely not be made completely in the first operation (errors!). They are done step by step (usually every second day by a maximum of 2 mm!).

Durchführung der **Ilizarov**-Ringfixateurmontage Mounting the **Ilizarov** ring fixing device

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In supinated pointed foot, we first correct the rear foot (less stretching of the joints makes the correction possible!) and then the frontal foot. On the frontal foot, we first correct supination, and then gradually the pointed foot position.

If a foot malposition must be corrected, it is always necessary to insert at least a frontal foot wire aside from the heel ring (Fig. 10a). In supination foot, the olive is placed medially on the frontal foot. The wire from the first metatarsal bone is brought directly into the fifth metatarsal bone, that is, angled in both levels; one aims at the protuberance, wherein the assistant must push the metatarsal bones 2-4 into the ventral direction while the wire is being drilled! One does not attempt fixation from the first to the fifth metatarsal bones = flat foot, but, insofar as possible, only the first and fifth bones. After some practice on the model, this will be successful, however it is usually the most difficult wire. Malpositions are not tolerated later and always cause problems!! (Fig. 10a)

The frontal foot wire is fastened on a sufficiently large half ring in order to prevent problems in swelling and allow a bandage. (Fig. 10a+10b)

The half ring is connected to the lower leg with various small parts (Fig. 10b+10c) so that gradual correction of malpositions is possible. In the case of contractile pointed foot, the correction should only be initiated during surgery – but not done completely (skin necroses, severe pain, emerging wires).

Even complex foot malpositions can be gradually corrected with this foot fixation using 3 wires with corresponding rings.

The completed lower leg assembly with foot inclusion is shown in Fig. 10b-10d.

Durchführung des Segmenttransports beim **Ilizarov**-Ringfixateurmontage
Performance of segment transport with the **Ilizarov** ring fixing device

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To perform segment transport, it is necessary to create a transport segment on the one hand and one or more drawing mechanisms on the other hand.

The following can be used as drawing mechanisms (Fig. 11a):

1. Transport ring

This is a further ring between the second and third fixation rings which is transported into the distal or proximal direction on threaded rods which are precisely parallel.

The transport segment is captured with two opposing olive K-wires; the wires “cut” through the soft tissues and leave so-called scar tracks (which are usually not very noticeable).

2. Crossed olive K-wires (Fig. 9a)

In the original Ilizarov method, olive K-wires or wires with ends which had loops or spindles were always used for transport. The wires are drilled through the tip of the transport segment at an angle and then driven through the soft tissues, brought out distally or proximally (depending on the pull direction) and fastened on threaded rods with slits (Fig. 11b), which in turn are fastened with 3 end pieces (all levels of freedom) with thread on a ring or other auxiliary structures.


The threaded rods with slits can then be gradually turned up without problems, thereby enabling transport.

3. Cerclage wires, which are anchored on the transport segment with screws.

To improve transport (olive K-wires are extremely rigid), we introduced cerclage wires which are fastened on the thigh with 2 and on the lower leg with one screw (see Fig. 11a).

The cerclage wires are intraoperatively brought out through the soft tissues, and then gradually moved using the same drawing mechanisms as in olive K-wires.

Bringing the cerclage wires through the soft tissues is simplified by bringing out redon skewers with redon from the defect through the soft tissues where the segment will be transported. Then one simply places the wire in the redon, making it possible to place them in the correct exit point without problems.

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Performance of segment transport with the **Ilizarov** ring fixing device

4. Rope system according to Weber

Weber [Weber 2003] introduced a rope pulley system by using the pulley ropes in place of cerclage wires (more flexible) on the one hand, and on the other hand, using redirecting rollers which make it possible to redirect the distal pull ropes into the proximal direction. The patient can perform the repeated daily turning much more easily into the proximal direction than into the distal direction.

The rope pulley system is performed through two angled drill holes with a diameter of 2 or 2.7 mm at the distal end of the transport segment in a figure eight from one side, so that it is not necessary to bare a large section of bone. When the transport is complete, the rope is simply pulled out (Fig. 11a).

5. Corticotomy

The so-called corticotomy produces the transport segment.

The corticotomy (bone separation while maintaining the endostal vessels) can be done in various ways. We are meanwhile performing corticotomies after carrying out several fan shaped drill holes on the bone with Ilisarow's special chisel, wherein we always attempt to maintain endostal circulation; naturally, this is not always possible.

The corticotomy must be performed completely. If an unstable bone connection remains on only one side, the segment transport is not successful.

When carrying out the segment transport, many special attributes which cannot be addressed here in further detail must be taken into account (see literature).

This contribution merely shows the fundamental implementation of ring assembly. Using the ring system in accordance with Ilizarow, it is possible to build many static, dynamic or combined assemblies. Finally, it must be emphasized that the forming distraction bone – due to its intramembranous formation – leads primarily to the tubular bone and does not require restructuring, which is a significant advantage as compared to restructuring after spongiosa transplantation.

Elemente Ilizarov Ringfixateur
elements ilizarov ring fixator

Halbring half ring				
	Ø Ø	Dicke thickness	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	60 mm	4 mm	1001060T	
	70 mm	4 mm	1001070T	
	80 mm	4 mm	1001080T	10.01.080
	100 mm	4 mm	1001100T	10.01.100
	110 mm	4 mm	1001110T	10.01.110
	120 mm	4 mm	1001120T	10.01.120
	130 mm	4 mm	1001130T	10.01.130
	140 mm	4 mm	1001140T	10.01.140
	150 mm	4 mm	1001150T	10.01.150
	160 mm	4 mm	1001160T	10.01.160
	180 mm	5 mm	1001180T	10.01.180
	200 mm	5 mm	1001200T	10.01.200
	220 mm	5 mm	1001220T	10.01.220
	240 mm	5 mm	1001240T	10.01.240

5/8 Ring 5/8 ring				
	Ø Ø	Dicke thickness	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	130 mm	4 mm	1003130T	10.03.130
	150 mm	4 mm	1003150T	10.03.150
	160 mm	4 mm	1003160T	10.03.160

Vollring full ring				
	Ø Ø	Dicke thickness	Art. Nr. Titan ref. no. titanium	Art. Nr. ref. no. steel
	Stahl			
	100 mm	4 mm	1002100T	10.02.100
	110 mm	4 mm	1002110T	10.02.110
	120 mm	4 mm	1002120T	10.02.120
	130 mm	4 mm	1002130T	10.02.130
	140 mm	4 mm	1002140T	10.02.140
	150 mm	4 mm	1002150T	10.02.150
	160 mm	4 mm	1002160T	10.02.160
	180 mm	5 mm	1002180T	10.02.180
	200 mm	5 mm	1002200T	10.02.200
	220 mm	5 mm	1002220T	10.02.220
	240 mm	5 mm	1002240T	10.02.240

Elemente Ilizarov Ringfixateur
elements ilizarov ring fixator

Bogenplatte curved plate			
	Biegung Curvature um um	Dicke thickness	Art. Nr. Titan ref. no. titani-
	90°	4 mm	1004090T
	120°	4 mm	1004120T

Halbring mit aufgebogenen Enden half ring with curved ends				
	Ø Ø	Dicke thickness	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	120 mm	4 mm	1005120T	10.05.120
	130 mm	4 mm	1005130T	10.05.130
	150 mm	4 mm	1005150T	10.05.150

Fußstütze foot ring			
	Länge length um um	Dicke thickness	Art. Nr. Titan ref. no. titani-
	260 mm	4 mm	1003260T
	300 mm	4 mm	1003300T

Gewindestange M6 threaded rod M6				
	Länge length	Anz. (VE) qty. (PU)	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	60 mm	5	1006060T	10.06.060
	80 mm	5	1006080T	10.06.080
	100 mm	5	1006100T	10.06.100
	120 mm	5	1006120T	10.06.120
	150 mm	5	1006150T	10.06.150
	200 mm	5	1006200T	10.06.200
	250 mm	5	1006250T	10.06.250
	300 mm	5	1006300T	10.06.300
	350 mm	5	1006350T	10.06.350
	400 mm	5	1006400T	10.06.400

Gewindestange M6 mit Schlitz threaded rod M6 slotted				
	Länge length	Anz. (VE) qty. (PU)	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	40 mm	5	1007040T	10.07.040
	60 mm	5	1007060T	10.07.060
	80 mm	5	1007080T	10.07.080
	100 mm	5	1007100T	10.07.100
	150 mm	5	1007150T	

Elemente Ilizarov Ringfixateur
elements ilizarov ring fixator

Verbindungsplatte kurz connection plate short				
	Länge length	Loch hole	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	37 mm	2	1010037T	
	49 mm	3	1010049T	10.10.049
	60 mm	4	1010060T	10.10.060
	71 mm	5	1010071T	10.10.071
	82 mm	6	1010082T	10.10.082
	93 mm	7	1010093T	10.10.093
	104 mm	8	1010104T	10.10.104
	115 mm	9	1010115T	10.10.115
	126 mm	10	1010126T	10.10.126

Verbindungsplatte lang connection plate long				
	Länge length	Loch hole	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	159 mm	14	1011159T	10.11.159
	236 mm	21	1011236T	10.11.236
	335 mm	30	1011335T	10.11.335

Verbindungsplatte mit Gewindeansatz connection plate with threaded end				
	Länge length	Loch hole	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	131 mm	5	1012131T	10.12.131
	171 mm	7	1012171T	10.12.171
	211 mm	9	1012211T	10.12.211
	251 mm	11	1012251T	10.12.251

Verwindungsplatte twisted plate				
	Länge length	Loch hole	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	46 mm	2	1013046T	10.13.046
	67 mm	3	1013067T	10.13.067
	88 mm	4	1013088T	10.13.088

Endteil mit Innengewinde endpiece with inner thread		
	Loch hole	Art. Nr. Titan ref. no. titanium
	1	1015001T
	2	1015002T
	3	1015003T
	4	1015004T

Endteil mit Gewindeansatz endpiece with threaded end		
	Loch hole	Art. Nr. Titan ref. no. titanium
	1	1016001T
	2	1016002T
	3	1016003T
	4	1016004T

Elemente Ilizarov Ringfixateur
elements ilizarov ring fixator

Gewindemuffe Sechskant M6 threaded sleeve hexagonal M6		
	Länge length	Art. Nr. Titan ref. no. titanium
	20 mm	1018020T
	30 mm	1018030T
	40 mm	1018040T
	60 mm	1018060T

Drahtfixationsklemme wire fixation clamp		
	Art. Nr. Titan ref. no. titanium	Art. Nr. Stahl ref. no. steel
	1022000T	10.22.000

Unterlegscheibe für Olivendraht washer for olive wire		
	Anz. (VE) qty. (PU)	Art. Nr. Titan ref. no. titanium
	5	1027007T

Segmenttransporteinheit segment transport unit	
	Art. Nr. Stahl ref. no. steel
	10.28.000

Titanseil titanium cable			
	Ø Ø	Länge length	Art. Nr. Titan ref. no. titanium
	1,15 mm	5 m	16115100T

Umlenkrolle für Segmenttransport pulley block for segment transport		
	Ø Ø	Art. Nr. Titan ref. no. titanium
	15 mm	1021015T

Elemente Ilizarov Ringfixateur

elements ilizarov ring fixator

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safety washer

Art. Nr. Stahl

ref. no. steel

10.21.006

Unterlegscheibe

washer

Dicke

thickness

Anz. (VE)

qty. (PU)

Art. Nr.

ref. no.

1 mm

10

1021001T

2 mm

10

1021002T

3 mm

(mit Schlitz)

(with slot)

10

1021003T

4 mm

10

1021004T

Kugelscheiben

rocker washer

Anz. (VE)

qty. (PU)

Art. Nr. Titan

ref. no. titanium

Art. Nr. Stahl

ref. no. steel

paarweise / pairwise

1021005T

10.21.012

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steinmann nail with middle flat thread

Ø Nagel

Ø nail

Ø Gewinde

Ø thread

Länge

length

Art.-Nr. Titan

ref.-no. titani-
um

4,5 mm

5,0 mm

175 mm

1036175T

4,5 mm

5,0 mm

200 mm

1036200T

5,0 mm

5,5 mm

175 mm

1035175T

5,0 mm

5,5 mm

200 mm

1035200T

5,0 mm

5,5 mm

250 mm

1035250T

5,0 mm

5,5 mm

300 mm

1035300T

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steinmann nail with middle sharp thread

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WIRE with bayonet point and stopper

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ABBOTT LC (1927)
The operative lengthening of the tibia and fibula. J Bone Joint Surg [Am] 9:128-152

ALLAN FG (1948)
Bone lengthening. J Bone Joint Surg [Br] 30:490

ANDERSON MV
(1952) Leg lengthening. J Bone Joint Surg [Br] 34:150

ARONSON J, HARP JH (1992)
Mechanical considerations in using tensioned wires in a transosseous external fixation system. Clin Orthop 280:23-29

ARONSON J, HARRISON B, BOYD CM, CANNON DJ, LUBANSKY HJ (1988)
Mechanical induction of osteogenesis: The importance of pin rigidity. J Pediatr Orthop 8: 396-401

ARONSON J, HARRISON BH, STEWART CL, HARP JH (1989)
The histology of distraction osteogenesis using different external fixators. Clin Orthop 241:106-116

BARR JS, OBER FR (1933)
Leg lengthening in adults. J Bone Joint Surg [Am] 15: 674-678

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DELLOYE C, DELEFORTRIE G, COUTELIER L, VINCENT A (1990)
Bone regenerate formation in cortical bone during distraction lengthening - An experimental study. Clin Orthop 250:34-42

DELPRETE C, GOLA MM (1993)
Mechanical performance of external fixators with wires for the treatment of bone fractures - Part I: Load-displacement behavior. J Biomech Eng 115:29-36

Literaturliste **Ilizarov** Ringfixateurmontage
Bibliography **Ilizarov** ring fixing

DELPRETE C, GOLA MM (1993)
Mechanical performance of external fixators with wires for the treatment of bone fractures - Part II: Wire tension and slippage. J Biomech Eng 115:37-42

GREEN SA (1991)
The Ilizarov method: Rancho technique. Orthop Clin North Am 22:677-688

ILIZAROV GA (1989)
The tension-stress effect on the genesis and growth of tissues. Part I: The influence of stability of fixation and soft-tissue preservation. Clin Orthop 238:249-281

ILIZAROV GA (1989)
The tension-stress effect on the genesis and growth of tissues: Part II: The influence of the rate and frequency of distraction. Clin Orthop 239:263-285

KAWAMURA B, HOSONO S, TAKAHASHI T, YANO T, KOBAYASHI Y, SHIBATA N, SHINODA Y (1968)
Limb lengthening by means of subcutaneous osteotomy. J Bone Joint Surg [Am] 50:851-878

KENWRIGHT J, WHITE SH (1993)
A historical review of limb lengthening and bone transport. Injury, 24 Suppl 2:54-61

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MONTICELLI G, SPINELLI R (1981)
Distraction epiphysiolysis as a method of limb lengthening. I.Experimantal study. Clin Orthop 154:254-261

PALEY D (1988)
Current techniques of limb lengthening. J Pediatr Orthop 8:73-92

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SHEARER JR, ROACH HI, PARSONS SW (1992)
Histology of a lengthened human tibia.
J Bone Joint Surg [Br] 74:39-44

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WHITE SH, KENWRIGHT J (1990)
The timing of distraction of an osteotomy. J Bone Joint Surg [Br] 72:356-361

Literaturliste **Ilizarov** Ringfixateurmontage
Bibliography **Ilizarov** ring fixing

WHITE SH, KENWRIGHT J (1991)
The importance of delay in distraction of osteotomies. Orthop Clin North Am 22:569-579
WOLFSON N, HEARN TC, THOMASON JJ, ARMSTRONG PF (1990)
Force and stiffness changes during Ilizarov leg lengthening. Clin Orthop 250:58-60

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technische Änderungen vorbehalten / für Druckfehler übernehmen wir keinerlei Haftung / copyright 12/2025 / Innovations Medical GmbH
Innovations Medical GmbH
subject to technical changes/ we assume no liability for printing errors / copyright 12/2025 / Innovations Medical GmbH